



F205 series rotator

- Rotary shaker for test tubes suitable for mixing different sized tubes via a “tumbling” motion and “rolling” motion
- Ideal to use in laboratories for preparation protocols
- Suitable for using with incubators, even refrigerated
- Adjustable angle rotation
- Digital speed control
- 0-90 min or ∞

| Technical data | | F205 |
|---------------------|------|--------------------|
| Speed | rpm | tra 2 e 40 |
| Timer | min | 0 ÷ 90 o ∞ |
| Display | | digital |
| Rotation control | | digital |
| Working temperature | °C | da +4 a +60 |
| Power | W | 80 |
| Dimension LxWxH | mm | 210 x270x240 |
| Power | V/Hz | 230/50 |
| Weight | kg | 3,2 |
| Item | | 603.0129.26 |

ACCESSORIES

Rotating plate tube holder

| Item | n. tubes | for ml |
|--------------------|----------|--------|
| 603.0129.41 | 48 | 1,5 |
| 603.0129.42 | 12 | 50 |

Rotating plate tube holder

| Item | n. tubes | Ø mm |
|--------------------|----------|------|
| 603.0129.43 | 20 | 9÷20 |

Rack tube holder for cell culture

| Item | n. tube | Ø mm |
|--------------------|---------|------|
| 603.0129.44 | 63 | 12 |
| 603.0129.45 | 63 | 16 |
| 603.0129.46 | 30 | 26 |



Roller

- Provides a rocking and rolling motion to test tubes
- Digital speed control
- 540 min digital timer
- Digital speed control with a resolution of 1 rpm
- Audible and visual alarm for motor overload
- Can be used with Bijoux, Universal containers, and bottles with diameters up to 135 mm

| Technical data | | F406 | F409 |
|-----------------------|------|--------------------|--------------------|
| Speed | rpm | 5-60 | 5-60 |
| Power | W | 50 | 50 |
| Timer | min | 540 | 540 |
| Amplitude excursion | mm | 16 | 16 |
| Rollers | n | 6 | 9 |
| Maximum load | kg | 10 | 10 |
| Internal dimension | mm | 340x30 | 340x30 |
| Total dimension LxWxH | mm | 565x240x110 | 565x360x110 |
| Environment | °C | da +4 a +60 | da +4 a +60 |
| Operating temperature | % | 80 | 80 |
| Power | V/Hz | 230/50 | 230/50 |
| Weight | kg | 5,1 | 6,9 |
| Item | | 603.0129.80 | 603.0129.82 |

Vortex mixer



Mix

- Ideal for mixing solutions in test tubes. It functions automatically by exerting pressure on the small stirring cup.
- Continuous operation can be selected by pushing the switch on the front panel. Stirring frequency can also be selected by the knob on the front panel.
- The Mixer is fitted with four small suction cups providing a firm grip on the work to.

| Technical data | | Mix |
|------------------|------|--------------------|
| Rotary motion | mm | 4+ |
| Variable speed | rpm | 100-3000 |
| Safety class | | 0 |
| Protection class | IP | 41 |
| Motor Power | W | 30 |
| Dimension LxWxH | mm | 145×175×140 |
| Power | V/Hz | 220/240 – 50/60 |
| Weight | Kg | 1,9 |
| Item | | 603.0129.02 |

ACCESSORIES



Clamp

Item

603.0129.06



Spring lock tube

Item

603.0129.10



Rod

Item

601.0125.01



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2 | 24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it | falcinstruments.it