OT 012

BENCH TOP STEAM STERILIZER



OT 012 bench top steam sterilizer is designed for the sterilization of unwrapped medical or dental instruments and materials used in laboratories with pressurized steam. By means of the programmable microprocessor control system, sterilization temperature and time could be easily programmed and there are digital displays for temperature and time. In case of electricity interruption, if the chamber temperature is in the limits of programmed alarm limits, the program continues to operate from the point it is interrupted. But, if the chamber temperature is out of the programmed alarm limits, the program is cancelled and the user is warned by an error message.

miras

The chamber pressure is observed from the manometer on the control panel. There is a safety valve against over pressure and the energy of the heater is cut in case of over pressure.

The operation of OT 012 is very simple and

there is no need for any water or drain connection.



TECHNICAL SPECIFICATIONS

OT 012
110°C / 140°C
I°C
± 1°C
I minute / 100 hours
Pressure switch and safety valve
230 V, 50 Hz
1500 W
12 lt.
1 / 4
Stainless steel
Epoxy-polyester painted stainless steel
200×330×200
400×500×400
38 / 42
560×630×440



NÜVE SANAYİ MALZEMELERİ **IMALAT VE TİCARET A.Ş.**

PRESSURE

Esenboğa Yolu, 22 km.
Akyurt 06287 ANKARA TURKEY
Tel: (90.312) 399 28 30-31-32 - 399 24 87
Fax: (90.312) 399 21 97
http://www.nuve.com.tr
e-mail: nuveexp@nuve.com.tr

ISO 9001 ISO 13485



